

DXP



DXP DOMINA

Effective reliable sterilization.

DXP Domina is designed in every detail to guarantee effective class B sterilisation, making the operator work in comfort and tranquillity and ensuring maximum reliability. The large copper chamber is characterised by its excellent loading capacity and the rapid and uniform heating; the Adaptive Heat system eliminates temperature fluctuations and allows even sensitive instruments while minimising the risk of damaging them. Thanks to the printer and optional Wi-Fi Kit, the data from each cycle is stored easily and accurately, for maximum peace of mind and safety of operators and patients.

Distributed by:



QUANTUM
BIOMEDICAL

Quantum Biomedical Services Sdn. Bhd.

HQ : Unit L3A-6, Level 3A, Wisma Kemajuan, No. 2 Jalan 19/1b, 46300 Petaling Jaya, Selangor Darul Ehsan.

Tel: 03-7629 6777 - Fax: 03-7629 6779 - Email: enquiries@qbms.com.my - Website: qbms.com.my

DXP DOMINA

STERILIZATION CYCLES

Cycle	Temperature	Sterilization	Drying	Vacuum	Type of Cycle	Maximum
B134	134° C	5'	10'	3	B	4 kg solid 1.5 kg porous
B121	121° C	20'	11'	3	B	4 kg solid 1.5 kg porous
S134	134° C	3'	5'	2	S	4 kg solid
CUSTOM	105 - 135° C	3' - 9'	5' - 12'	2 - 3 - 4		4 kg solid 1.5 kg porous
HELIX TEST	134° C	3' 30"	10'	3	TEST	
B&D TEST	134° C	3' 30"	10'	3	TEST	
VACUUM TEST	< 35° C			3	TEST	

TECHNICAL CHARACTERISTICS

	DXP DOMINA - 18 L
OVERALL DIMENSIONS	447 x 438 x 625 mm (L x H x D)
CHAMBER DIMENSIONS	240 - 384 mm (Ø x D)
CHAMBER VOLUME	17.5 litres
WEIGHT (empty)	45 kg
POWER SUPPLY VOLTAGE	230 V 50/60 Hz
CONSUMPTION (power)	900 W - 1920 W
OPERATING TEMPERATURE	+10 - +40°C
OPERATING ALTITUDE (leaving the loading cap free)	between 0 and 2000 metres above sea level
OPERATING RANGE	Indoor use Relative humidity 85% Max. mains voltage variation ±10% Installation category II Pollution level 2
VOLUME OF WATER CONSUMED	0.6 litres
TANK CAPACITY (clean water and used water)	4 litres each
AUTOMATIC SWITCH-OFF	After 30' on standby at the end of the cycle



Distributed by:



**QUANTUM
BIOMEDICAL**

Quantum Biomedical Services Sdn. Bhd.

HQ : Unit L3A-6, Level 3A, Wisma Kemajuan, No. 2 Jalan 19/1b, 46300 Petaling Jaya, Selangor Darul Ehsan.

Tel: 03-7629 6777 - Fax: 03-7629 6779 - Email: enquiries@qbms.com.my - Website: qbms.com.my



DXP - NSK Dental Italy S.r.l. - 36016 Thiene (VI) - Italy - Tel. +39 0445 820070 - www.dxp-sterilization.com